

Single Channel Pipettes – Finnpipette® Fixed Volume

A fixed volume pipette is ideal for repeat pipetting of a single volume. Fixed volume Finnpipettes offer the same features and advantages as other Finnpipette Digital pipettes.

Finnpipette Fixed Volume pipettes are available in the standard sizes listed below. Non-standard volumes are available on request, at no extra charge.

Finnpipette Fixed Volume						
Cat. no.	Volume μ l	Accuracy		Precision		Fintip
		μ l	%	s.d.* μ l	CV %*	
4501 000	1 μ l	± 0.04	± 4.0	0.04	4.0	10, 20, 50
4501 010	2 μ l	± 0.05	± 2.5	0.04	2.0	10, 20, 50
4501 020	5 μ l	± 0.07	± 1.4	0.07	1.4	10, 20, 50
4501 030	10 μ l	± 0.09	± 0.9	0.08	0.8	10, 20, 50
4501 040	20 μ l	± 0.12	± 0.6	0.10	0.5	250 Universal, 300, 200 Ext
4501 050	25 μ l	± 0.15	± 0.6	0.13	0.5	250 Universal, 300, 200 Ext
4501 060	50 μ l	± 0.3	± 0.6	0.2	0.4	250 Universal, 300, 200 Ext
4501 070	100 μ l	± 0.4	± 0.4	0.3	0.3	250 Universal, 300, 200 Ext
4501 080	200 μ l	± 0.8	± 0.4	0.6	0.3	250 Universal, 300, 200 Ext
4501 090	250 μ l	± 1.0	± 0.4	0.8	0.3	1000, 1000 Ext
4501 100	500 μ l	± 1.5	± 0.3	1.5	0.3	1000, 1000 Ext
4501 110	1000 μ l	± 3.0	± 0.3	3.0	0.3	1000, 1000 Ext
4501 120	2000 μ l	± 6.0	± 0.3	4.0	0.2	5 ml
4501 130	3000 μ l	± 9.0	± 0.3	6.0	0.2	5 ml
4501 140	5000 μ l	± 15.0	± 0.3	10.0	0.2	5 ml
4501 150	10 ml	± 30.0	± 0.3	20.0	0.2	10 ml

Please note: 0.2 μ l to 10 ml models are available on request.
* s.d. = Standard Deviation, CV = Coefficient of Variation

- Single volume repeating
- Non-standard volumes available
- Soft-touch tip ejection
- Fully autoclavable
- Ergonomic design



13

[Volver a productos >>](#)