

# HAROLDAN<sup>®</sup> Plus



**Safety lancets for capillary blood sampling**

### Features and benefits

- Needle retracts automatically and locks out for extra safety, protecting patient and health professionals from needlestick injuries and prevent cross-infection.
- Ergonomic and easy to handle.
- HAEMOLANCE Plus saves time. It is quick and simple to use and requires little instruction.
- HAEMOLANCE Plus can not be reloaded or reused.
- Simple and safe for both patient and user!

### How to use HAEMOLANCE Plus



1. To remove the protective cap turn and twist in one movement (a full turn).



2. Place lancet against side of fingertip and press button firmly. Discard used lancet into a sharps container.

	Name	Needle Size	Puncture Depth	Product No.	Pcs/Box	Transport and ordering pack.
	Paediatric	Blade	1,2 mm	7564	150 pcs	24 boxes
	Max Flow	Blade	1,6 mm	7591	150 pcs	24 boxes
	High Flow	≈18G	1,8 mm	7586 7585	150 pcs 50 pcs	24 boxes 54 boxes
	Normal Flow	≈21G	1,8 mm	7568 7567	150 pcs 50 pcs	24 boxes 54 boxes
	Low Flow	≈25G	1,4 mm	7575 7574	150 pcs 50 pcs	24 boxes 54 boxes
	Micro Flow	28G	1,6 mm	7561 7560	150 pcs 50 pcs	24 boxes 54 boxes

